alfresa

FIT Hemoglobin AA01

Product list

Catalogue No.	Product name
912240	FIT Hemoglobin AA01

Manufactured by:

Alfresa Pharma Corporation

18 Taiheidai Shouou-Cho, Katsuta-gun, Okayama, 709-4321, Japan

European Authorized Representative

Emergo Europe

FOR *IN VITRO* DIAGNOSTIC USE ONLY

1.0 INTENDED USE

1.1 **FIT Hemoglobin AA01** is a reagent to quantitatively determine human hemoglobin concentrations in feces using the **Discrete Clinical Chemistry Analyzer AA01**, which employs a colloidal gold immune colorimetric and turbidimetric method. This reagent is intended for colorectal cancer screening. FIT stands for fecal immunochemical test.

2.0 PRINCIPLE OF THE EXAMINATION METHOD

2.1 Summary and explanation of the test

Fecal immunochemical test (FIT) is used to diagnose hemorrhagic gastrointestinal diseases. FIT is particularly useful for colorectal cancer screening. **FIT Hemoglobin AA01** is a kit to measure human hemoglobin concentrations in feces using an immunochemical method combined with a colloidal gold colorimetric and turbidimetric method. This colloidal gold immune colorimetric and turbidimetric method is intended to measure an optical color change due to agglutination between colloidal gold-conjugated rabbit anti-human hemoglobin polyclonal antibodies and fecal human hemoglobin. This test is highly specific and sensitive.

2.2 Principle of the test

The reaction of colloidal gold-conjugated anti-human hemoglobin polyclonal antibodies with human hemoglobin in feces produces a color change due to agglutination of colloidal gold particles through the antigen-antibody reaction. Human hemoglobin concentration in feces is determined by measuring the color change over time.

3.0 COMPONENTS

3.1 FIT Hemoglobin AA01

A kit contains two sets of Reagent 1 and Reagent 2 bottles. Each bottle is for 125 tests. Total 250 tests can be assayed per kit.

3.1.1 Reagent 1 REAG1

FIT Hemoglobin AA01 contains 2 bottles of 11mL of Reagent 1. Each bottle contains:

MES Buffer	150 mmol/L
Sodium azide	< 0.1 %

Reagent 1 contains less than 0.1 % sodium azide. Attention, see **7.0 WARNINGS AND PRECAUTIONS**. Material Safety Data Sheet is available upon request by a professional user.

3.1.2 Reagent 2 REAG2

FIT Hemoglobin AA01 contains 2 bottles of 5 mL of Reagent 2. Each bottle contains:

TES Buffer	4.8 mmol/L
Bovine serum albumin	0.03 %
Sodium azide	< 0.1 %
Colloidal gold-conjugated anti-human hemoglobin polyclonal antibodies	267 µL/mL

Reagent 2 contains less than 0.1 % sodium azide. Attention, see **7.0 WARNINGS AND PRECAUTIONS**. Material Safety Data Sheet is available upon request by a professional user.

4.0 ADDITIONAL REQUIRED EQUIPMENT

4.1 Analyzer

Discrete Clinical Chemistry Analyzer AA01

4.2 Specimen Collection Container

Specimen Collection Container A

4.3 AA01 dedicated reagent

FIT Hemoglobin NS-Prime Calibrator

FIT NS-Prime Control

FIT AA01 Specimen Diluent

Wash Solution A

5.0 REAGENT PREPARATION

Mix the reagent 2 thoroughly before each use, because colloidal gold particles may precipitate during storage. Mix it by slowly inverting to avoid creating bubbles. If bubbles occur, the instrument sensor will read an incorrect volume and the bottle will be deleted from the system. Pipette out any bubbles down to the level of the liquid before placing the bottle on board.

6.0 STORAGE AND SHELF LIFE AFTER FIRST OPENING

6.1 Storage

Store at 2–8 °C. Do not freeze.

6.2 Storage and shelf life after first opening

Once opened, use the reagent within one month and store at 2–8 °C.

7.0 WARNINGS AND PRECAUTIONS

7.1 General precautions

For *In vitro* diagnostic use

Only experienced laboratory personnel should use this; the test should be used in a manner consistent with Good Laboratory Practice.

7.2 Safety precautions

- Do not pipet by mouth.
- Some reagents contain less than 0.1 % sodium azide. Upon exposure to the eye or skin or accidental ingestion, take emergency measures such as washing with plenty of water. Consult a doctor if necessary.
- Do not smoke, eat, or apply cosmetics in areas where patients' specimens or kit reagents are handled.
- Cuts, abrasions, and other skin lesions should be properly protected with an appropriate waterproof dressing.
- Take care to avoid self-inoculation, splashing to the mucous membrane, or generation of aerosols.
- Wear laboratory gloves while handling patients' specimens or disposing of solid or liquid wastes.
- Cautions upon disposal
 - Reagent 2 of **FIT Hemoglobin AA01** contains 0.22 g/L ethylenediaminetetraacetic acid copper (II) disodium (30 mg/L as copper). Upon disposal, comply with relevant legal provisions.
 - Some reagents contain less than 0.1 % sodium azide. Sodium azide can react with copper and lead plumbing to form explosive metal azides. Regulations currently in use regarding dangerous waste elimination must be followed. If disposed in the sink, rinse with plenty of water.
 - Upon disposal of reagents or other materials, comply with relevant legal provisions.
- Some reagents contain bovine serum albumin free from known infectious agents. However they should be considered potentially infectious and handled with care to avoid infection.
- All human specimens should be considered potentially infectious. Handle all specimens as if capable of transmitting HBV, HCV, HIV, or other microbes. Decontaminate and dispose of specimens and all potentially contaminated materials as if they contain infectious agents.

7.3 Limitations

- Do not use reagent containers for purposes other than this test.
- Do not damage or stain the bar codes on the labels of each container.
- Do not replenish or mix reagents. Also, do not mix reagents of different bottles even if they have the same lot number.
- Do not use combinations of different lot numbers within the products.
- Do not recycle the bottle. It may be infectious.
- Use of reagents, disposables, or spare parts other than those supplied by the authorized distributor may produce incorrect results.

8.0 PRIMARY SAMPLE COLLECTION, HANDLING, AND STORAGE

Use human fecal specimens for measurement.

- Scrape the surface of fecal specimen with the stick of **Specimen Collection Container A** thoroughly to prepare a sample.
- Specimens whose concentrations exceed the upper limit of the calibration curve should be diluted with **FIT AA01 Specimen Diluent** and retested.

8.1 Specimen collection

Sample feces onto the gutters of the collector stick by scraping several surfaces of fecal specimen. Place the stick into the collector body once only and fasten the stick tightly. See the instruction manual for **Specimen Collection Container A**.

8.2 Specimen storage

The collected sample should be refrigerated between (2–8 °C) until the sample can be tested. Once the sample is received for testing store between 2–8 °C.

8.3 Hemoglobin stability after feces sampling

The stability of human hemoglobin in the buffer of **Specimen Collection Container A** was examined. After storage of four different concentrations of human hemoglobin for 7, 14, and 33 days at –40, 7, 25, and 37 °C, the residual ratios of human hemoglobin were as follows:

	–40 °C, for 33 days	7 °C, for 33 days	25 °C, for 14 days	37 °C, for 7 days
50 ng/mL	> 90 %	> 90 %	> 90 %	> 90 %
169 ng/mL	> 90 %	> 90 %	> 90 %	> 70 %
266 ng/mL	> 90 %	> 90 %	> 70 %	> 50 %*
381 ng/mL	> 90 %	> 90 %	> 50 %*	> 30 %*

*Hb residual ratios decreased to 37–55 %, but the results were still positive (cut-off: 100 ng/mL). This result is provided for reference only. The stability varies depending on the specimen.

9.0 EXAMINATION PROCEDURE

9.1 Preparation of specimens

- Sample feces using **Specimen Collection Container A**. See the instruction manual for the container.
- Shake the container sufficiently to dissolve the feces from the gutters of the stick.
- Leave the container at room temperature for 1 hour. When measuring on the next day, refrigerate at 2–8 °C. Depending on the specimen, dissolution of hemoglobin from feces may be insufficient in 1 hour.

9.2 Assay

- Prepared specimen solution

Reagent 1	9 µL
Reagent 2	75 µL
	30 µL

Mix the reaction liquid and incubate at 37 °C. Formula to calculate the change of Absorbance (Abs) is (Am1–As1)–(Am2–As2).
Am1 : Abs of main wavelength at measurement point 1
As1 : Abs of sub wavelength at measurement point 1
Am2 : Abs of main wavelength at measurement point 2
As2 : Abs of sub wavelength at measurement point 2

main wavelength: 525 nm
sub wavelength: 660 nm
measurement point 1: 0 minutes
measurement point 2: 6.7 minutes

- Construct a calibration curve using **FIT Hemoglobin NS-Prime Calibrator** in place of the specimen solution. Read the human hemoglobin concentration in the specimen solution from the calibration curve.

All routine maintenance procedures defined in the instrument manual for the **Discrete Clinical Chemistry Analyzer AA01** should be performed to obtain optimal performance.

10.0 CONTROL PROCEDURE

It is recommended that each laboratory should use quality control materials routinely such as **FIT NS-Prime Control**, and should establish its own control ranges. Multilevel controls should be tested for each run. The human hemoglobin values obtained for the quality control materials should not fall repeatedly outside the control ranges established in each laboratory. If these control values fall repeatedly outside of the established control ranges, then proper instrument performance should be verified or perform recalibration.

11.0 CALCULATION OF EXAMINATION RESULTS

Specimen values are calculated by the **Discrete Clinical Chemistry Analyzer AA01**.

12.0 INTERPRETATION OF RESULTS

The cut-off value widely used in Japan for mass-screening is less than 100 ng/mL (equivalent to 20 µg of human hemoglobin per 1 g of feces). Colorectal cancer is suspected if human hemoglobin is detected above that level.

<In-house data>

Note : This value is indicative only and may differ from other published values because of differences in methods and in the population being studied. It is recommended that each laboratory establish its own cut-off value.

13.0 PERFORMANCE CHARACTERISTICS

The following performance data was obtained using the **Discrete Clinical Chemistry Analyzer AA01**.

13.1 Analytical performance characteristics

13.1.1 Precision

Precision: CV ≤ 15 %

13.1.2 Interference

Some studies have been conducted to determine the levels of interference. The following substances do not affect the assay at the following concentrations:

Conjugated and unconjugated bilirubin: 0–20 µg/mL
Ascorbic acid: 0–0.1 mg/mL
Glucose: 0–0.1 mg/mL
Bovine serum albumin: 0–0.1 mg/mL
Peroxidase: 0–0.1 mg/mL
Barium sulfate: 0–1.0 mg/mL
The cross-reactivities of the reagents with hemoglobin in animal species other than humans (bovine, swine, equine, ovine, caprine, and leporine) did not exist at concentrations of 0–2.0 µg/mL.

13.1.3 Correlation

FIT Hemoglobin AA01 assay was compared with **FIT Hemoglobin NS-Prime** on 100 specimens. A linear regression analysis was performed and the following results were obtained:

x : **FIT Hemoglobin NS-Prime**
y : **FIT Hemoglobin AA01**
r = 0.995
y = 0.98x – 4.1

13.2 Diagnostic performance characteristics

13.2.1 Analytical sensitivity

The difference in the amounts of change in absorbance between hemoglobin concentrations of 0 and 100 ng/mL is 0.05 or above.

13.2.2 Analytical accuracy

Within 100±15 % of known concentration in the measurement of a control specimen of known concentration.

13.3 Measuring interval

13.3.1 Assay range

- Upper limit : 1,200 ng/mL (equivalent to 240 µg human hemoglobin per 1 g of feces)
- Lower limit
Limit of quantitation : 50 ng/mL (equivalent to 10 µg human hemoglobin per 1 g of feces)
Limit of detection : 20 ng/mL (equivalent to 4 µg human hemoglobin per 1 g of feces)

14.0 BIOLOGICAL REFERENCE INTERVALS

Less than 100 ng/mL (20 µg /g feces)

Note : This value is indicative only and may differ from other published values because of differences in methods and in the population being studied. It is recommended that each laboratory establish its own local reference ranges.

15.0 LIMITATION OF THE EXAMINATION PROCEDURE

- Cuts, abrasions, and other skin lesions should be properly protected with an appropriate waterproof dressing, because blood contamination may influence measurements.
- Appearance changes, such as cloudiness and aggregation, in any of the reagents indicate the possibility of deterioration. Call your local dealer for advice.
- Mix the reagent 2 thoroughly before each use, because colloidal gold particles may precipitate during storage. Mix it by slowly inverting to avoid creating bubbles. If bubbles occur, the instrument sensor will read an incorrect volume and the bottle will be deleted from the system. Pipette out any bubbles down to the level of the liquid.

4. As with all assays, the results of this test can be influenced by factors present in some patients' specimens.

- For diagnostic purposes, the results obtained from this assay should always be used in combination with a clinical examination, patient medical history, and other findings.
- Procedural directions must be followed exactly because any modification of the procedure may change the results.
- Use of reagents, disposables, or spare parts other than those supplied by the authorized distributor may produce incorrect results.
- Store the reagents according to the storage methods. Do not use them after the expiration date.
- Use fresh feces.
- This test should not be used to analyze specimens taken from a patient who is menstruating or who has hemorrhoids.
- The test has not been validated for testing of patients with hemoglobinopathies.

16.0 LITERATURE REFERENCES

- Yasuhiro Oono et al., A retrospective study of immunochemical fecal occult blood testing for colorectal cancer detection, *Clinica Chimica Acta*, 411, p802–805, 2010.
- Edward Randell et al., Evaluation of Hemo Techt NS-Plus system for use in a province-wide colorectal cancer screening program, *Clinical Biochemistry*, 46(4–5), p365–368, 2013.
- Setsuko Kato et al., Basic Evaluation of Hemotect NS-Prime, an Automated Immunochemical Analyzer for Fecal Occult Blood Testing, *J Clin Lab Inst Reag*, 37(3), p371–377, 2014.
- Wael L.L.Demian et al., Evaluation of the analytical performance of the novel NS-Prime system and examination of temperature stability of fecal transferrin compared with fecal hemoglobin as biomarkers in a colon cancer screening program, *Practical Laboratory Medicine*, 2, p29–36, 2015.
- Ari Ahn, M.D. et al., Performance Evaluation of Two Automated Quantitative Fecal Occult Blood Tests, *Lab Med Online*, 6(4), p233–239, 2016.
- Taira Inagaki et al., Evaluation of the new automatic immunochemistry fecal occult blood analyzer "Hemo Techt NS-Prime", *Japanese Journal of Medical Technology*, 65(2), p222–228, 2016.
- Takenobu Shimada et al., How should an adequate recall rate in colorectal cancer screening be evaluated and dealt with? , *Journal of Gastrointestinal Cancer Screening*, 55(1), p31–44, 2017.

17.0 SYMBOLS USED IN PRODUCT INSERTS AND ON LABELS

Symbols	Meanings of the symbols
	Expiry date (Used by…)
LOT	Lot number
REF	Catalogue code
	Manufactured by
EC REP	Authorized EC representative
	Number of tests
IVD	<i>In vitro</i> diagnostic medical device (<i>In vitro</i> diagnostic use)
	Temperature limitation (store at…)
	See instruction for use
	Infectious substances: Handle with care
REAG 1	Reagent 1
REAG 2	Reagent 2

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