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

17.0 LITERATURVERZEICHNIS

- Yasuhiro Oono et al., A retrospective study of immunochemical fecal occult blood testing for colorectal cancer detection, Clinica Chimica Acta, 411, p802–805, 2010.
- Edward Randell et al., Evaluation of Hemo Techt NS-Plus system for use in a province-wide colorectal cancer screening program, Clinical Biochemistry, 46(4-5), p365–368, 2013.
- Setsuko Kato et al., Basic Evaluation of Hemotect NS-Prime, an Automated Immunochemical Analyzer for Fecal Occult Blood Testing, J Clin Lab Inst Reag, 37(3), p371–377, 2014.
- Wael L.L.Demian et al., Evaluation of the analytical performance of the novel NS-Prime system and examination of temperature stability of fecal transferrin compared with fecal hemoglobin as biomarkers in a colon cancer screening program, Practical Laboratory Medicine, 2, p29–36, 2015.
- Ari Ahn, M.D. et al., Performance Evaluation of Two Automated Quantitative Fecal Occult Blood Tests, Lab Med Online, 6(4), p233–239, 2016.
- Taira Inagaki et al., Evaluation of the new automatic immunochemistry fecal occult blood analyzer "Hemo Techt NS-Prime", Japanese Journal of Medical Technology, 65(2), p222–228, 2016.
- Takenobu Shimada et al., How should an adequate recall rate in colorectal cancer screening be evaluated and dealt with? , Journal of Gastrointestinal Cancer Screening, 55(1), p31–44, 2017

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