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17.0 RIFERIMENTI BIBLIOGRAFICI

1. Yasuhiro Oono et al., A retrospective study of immunochemical fecal occult blood testing for colorectal cancer detection, Clinica Chimica Acta, 411, p802–805, 2010.
2. Edward Randell et al., Evaluation of Hemo Techt NS-Plus system for use in a province-wide colorectal cancer screening program, Clinical Biochemistry, 46(4-5), p365–368, 2013.
3. Setsuko Kato et al., Basic Evaluation of Hemotect NS-Prime, an Automated Immunochemical Analyzer for Fecal Occult Blood Testing, J Clin Lab Inst Reag, 37(3), p371–377, 2014.
4. Wael L.L.Demian et al., Evaluation of the analytical performance of the novel NS-Prime system and examination of temperature stability of fecal transferrin compared with fecal hemoglobin as biomarkers in a colon cancer screening program, Practical Laboratory Medicine, 2, p29–36, 2015.
5. Ari Ahn, M.D. et al., Performance Evaluation of Two Automated Quantitative Fecal Occult Blood Tests, Lab Med Online, 6(4), p233–239, 2016.
6. Taira Inagaki et al., Evaluation of the new automatic immunochemistry fecal occult blood analyzer "Hemo Techt NS-Prime", Japanese Journal of Medical Technology, 65(2), p222–228, 2016.
7. Takenobu Shimada et al., How should an adequate recall rate in colorectal cancer screening be evaluated and dealt with? , Journal of Gastrointestinal Cancer Screening, 55(1), p31–44, 2017

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