



## NESCAUTO Cp Auto Test

Catalogue No.	Product name
912681	NESCAUTO Cp Auto NS-Prime
912684	NESCAUTO Cp Auto Calibrator
912683	NESCAUTO Cp Auto Control

Note: NESCAUTO Cp Auto NS-Prime can only be used with Discrete Clinical Chemistry Analyzer NS-Prime. NESCAUTO Cp Auto Calibrator and NESCAUTO Cp Auto Control can be used with both Discrete Clinical Chemistry Analyzer NS-Prime and AA01.

#### Manufactured by:

Alfresa Pharma Corporation  
18 Taiheidai Shouou-Cho, Katsuta-gun, Okayama, 709-4321, Japan

#### European Authorized Representative

Emergo Europe  
Prinsessegracht 20 2514 AP The Hague The Netherlands

### Refer to the following URL for the package inserts in languages other than English:

<http://www.alfresa-pharma.co.jp/english/works/packageinserts.html>



### FOR *IN VITRO* DIAGNOSTIC USE ONLY

#### 1.0 INTENDED USE

**1.1 NESCAUTO Cp Auto NS-Prime** is a reagent to quantitatively determine human calprotectin concentrations in feces using **Discrete Clinical Chemistry Analyzer NS-Prime**, which employs a colloidal gold immune colorimetric and turbidimetric method. Fecal calprotectin levels reflect the degree of inflammation in the intestinal mucosa and are known to correlate with histologic and endoscopic evaluations.<sup>1,2)</sup> This reagent is intended as an aid in the assessment of intestinal mucosal inflammation in inflammatory bowel disease (IBD) patients, and as an aid in the differentiation of IBD from irritable bowel syndrome (IBS).

**1.2 NESCAUTO Cp Auto Calibrator** is a calibrator to construct a calibration curve on **Discrete Clinical Chemistry Analyzers NS-Prime** and **AA01**. Specimen calprotectin concentration is determined using the calibration curve.

**1.3 NESCAUTO Cp Auto Control** is a daily liquid control for routine quality assurance, which is used with **Discrete Clinical Chemistry Analyzers NS-Prime** and **AA01**. It should be run every day before specimen assay to confirm that its value is within the acceptable range.

#### 2.0 PRINCIPLE OF THE EXAMINATION METHOD

The calprotectin test is particularly useful for assessment of IBD patients. **NESCAUTO Cp Auto Test** is a kit to measure human calprotectin concentrations in feces using an immunochemical method combined with a colloidal gold colorimetric and turbidimetric method. This colloidal gold immune colorimetric and turbidimetric method is intended to measure the optical color change that takes place due to agglutination between colloidal gold-conjugated mouse anti-human calprotectin monoclonal antibodies and human calprotectin in feces. The color change is ascertained by measuring optical absorption using **Discrete Clinical Chemistry Analyzer NS-Prime** or **AA01**. Human calprotectin concentration in feces is calculated with the optical absorption change of specimen and the calibration curve generated from the measured calibrator values.<sup>3)</sup>

#### 3.0 TRACEABILITY OF VALUES ASSIGNED TO CALIBRATORS AND TRUENESS-CONTROL MATERIALS

Recombinant human calprotectin is a reference material to calibrate calibrators and controls. The reference value for recombinant human calprotectin is determined using the UV method.

#### 4.0 COMPONENTS

Caution: All components contain less than 0.1% sodium azide. Attention, see **8.0 WARNINGS AND PRECAUTIONS**. Safety Data Sheets are available upon request by professional users.

##### 4.1 NESCAUTO Cp Auto NS-Prime

Reagent 1 and Reagent 2 bottles are combined and packaged into two sets. One set provides 200 tests. A total of 400 tests can be conducted.

##### 4.1.1 Reagent 1 REAG1

Each bottle contains 31 mL solution. It contains:

HEPES Buffer	100 mmol/L
Sodium azide	< 0.1%

##### 4.1.2 Reagent 2 REAG2

Each bottle contains 11 mL solution. It contains:

TES Buffer	4.4 mmol/L
Bovine serum albumin	0.07%
Sodium azide	< 0.1%
Colloidal gold-conjugated mouse anti-human calprotectin monoclonal antibodies	333 $\mu$ L/mL

##### 4.2 NESCAUTO Cp Auto Calibrator

**NESCAUTO Cp Auto Calibrator** consists of four Calibrator bottles and one Calibrator Diluent bottle. The item contains the materials required for calibration to be carried out four times.

##### 4.2.1 Calibrator

Each bottle contains 1 mL solution. It contains:

MES Buffer	300 mmol/L
Sodium chloride	1.1%
Bovine serum albumin	0.15%
Sodium azide	< 0.1%
Human calprotectin	1200-1350 $\mu$ g/g (The assigned value is shown on the item label)

##### 4.2.2 Calibrator Diluent

Calibrator Diluent contains 12 mL solution. It contains:

MES Buffer	120 mmol/L
Sodium chloride	1.1%
Bovine serum albumin	0.15%
Sodium azide	< 0.1%

##### 4.3 NESCAUTO Cp Auto Control

**NESCAUTO Cp Auto Control** consists of Control L and Control H, and four of each are included.

##### 4.3.1 Control L

Each bottle contains 3 mL solution. It contains:

MES Buffer	300 mmol/L
Sodium chloride	1.1%
Bovine serum albumin	0.15%
Sodium azide	< 0.1%
Human calprotectin	120-200 $\mu$ g/g (The actual mean value and acceptable range are shown in the assigned values sheet.)

##### 4.3.2 Control H

Each bottle contains 3 mL solution. It contains:

MES Buffer	300 mmol/L
Sodium chloride	1.1%
Bovine serum albumin	0.15%
Sodium azide	< 0.1%
Human calprotectin	450-600 $\mu$ g/g (The actual mean value and acceptable range are shown in assigned values sheet.)

#### 5.0 ADDITIONAL REQUIRED EQUIPMENT

##### 5.1 Analyzer

**Discrete Clinical Chemistry Analyzer NS-Prime**

**Discrete Clinical Chemistry Analyzer AA01**

Note: **NESCAUTO Cp Auto NS-Prime** can only be used with **Discrete Clinical Chemistry Analyzer NS-Prime**, but **NESCAUTO Cp Auto Calibrator** and **NESCAUTO Cp Auto Control** can be used with both **Discrete Clinical Chemistry Analyzer NS-Prime** and **AA01**.

##### 5.2 Specimen Collection

**Specimen Collection Container A**

##### 5.3 Reagent

**NESCAUTO Cp Auto AA01**

Note: **NESCAUTO Cp Auto AA01** is used with **Discrete Clinical Chemistry Analyzer AA01**. It cannot be used with **Discrete Clinical Chemistry Analyzer NS-Prime**.

##### 5.4 Specimen Diluent

**FIT NS-Prime Specimen Diluent**

**FIT AA01 Specimen Diluent**

Note: **FIT NS-Prime Specimen Diluent** can only be used with **Discrete Clinical Chemistry Analyzer NS-Prime** and **FIT AA01 Specimen Diluent** can only be used with **Discrete Clinical Chemistry Analyzer AA01**.

##### 5.5 Wash Solution

**Wash Solution A**

##### 5.6 Others

Sample cup  
Plastic Pipette

#### 6.0 REAGENT PREPARATION

##### 6.1 NESCAUTO Cp Auto NS-Prime

As colloidal gold particles may precipitate during storage, mix the Reagent 1 and Reagent 2 combined bottle thoroughly before each use. Mix it by slowly inverting to avoid creating bubbles. If bubbles form, the instrument sensor cannot detect the liquid surface correctly and misunderstands as if the reagent volume increased (reagent from another bottle was added). As the result, the bottle will be deleted from the system. Pipette out any bubbles down to the level of the liquid before loading the bottle.

##### 6.2 NESCAUTO Cp Auto Calibrator

Mix Calibrator by slowly inverting before use. Dispense at least 1 mL of Calibrator Diluent into a sample cup for constructing a calibration curve with **Discrete Clinical Chemistry Analyzers NS-Prime** or **AA01**.

##### 6.3 NESCAUTO Cp Auto Control

Mix Control by slowly inverting before each use.

#### 7.0 STORAGE AND SHELF LIFE AFTER FIRST OPENING

##### 7.1 Storage and shelf life

Do not freeze reagents. The shelf life of unopened reagents at 2–8°C is 12 months: see the expiry date on the package label.

##### 7.2 Storage and shelf life after first opening

##### 7.2.1 NESCAUTO Cp Auto NS-Prime

Once opened, store the reagent at 2–8°C and use within 8 weeks.

##### 7.2.2 NESCAUTO Cp Auto Calibrator

The calibrator must be used on the same day that its bottle is opened.

##### 7.2.3 NESCAUTO Cp Auto Control

Once opened, store the control at 2–8°C and use within 3 weeks.

#### 8.0 WARNINGS AND PRECAUTIONS

##### 8.1 General precautions

For *in vitro* diagnostic use

Procedures should only be undertaken by experienced laboratory personnel; tests should be conducted in a manner consistent with Good Clinical Laboratory Practice.

##### 8.2 Safety precautions

- Do not pipet by mouth.
- Reagents contain less than 0.1% sodium azide. In case of accidental ingestion or exposure to eyes or skin, take emergency measures such as washing with plenty of water. Consult a doctor if necessary.
- Do not smoke, eat, or apply cosmetics in areas where patients' specimens or kit reagents are handled.
- Cuts, abrasions, and other skin lesions should be properly protected with an appropriate waterproof dressing.
- Take care to avoid self-inoculation, splashing to the mucous membranes, or generation of aerosols.
- Wear laboratory gloves while handling patients' specimens or disposing of solid or liquid waste.
- Cautions upon disposal
  - Reagent 2 of **NESCAUTO Cp Auto NS-Prime** contains 0.20 g/L ethylenediaminetetraacetic acid copper (II) disodium (27 mg/L as copper). Upon disposal, comply with relevant legal provisions.
  - Reagents contain less than 0.1% sodium azide. Sodium azide can react with copper and lead plumbing to form explosive metal azides. Regulations currently in use regarding dangerous waste elimination must be followed. If disposed of in the sink, rinse with plenty of water.
  - Upon disposal of reagents or other materials, comply with relevant legal provisions.
- Some reagents contain bovine serum albumin, free from known infectious agents. However, these reagents should be considered potentially infectious, and handled with care to avoid infection.
- All human specimens should be considered potentially infectious. Decontaminate and dispose of specimens and all potentially contaminated materials as if they contain infectious agents.
- NESCAUTO Cp Auto Calibrator** and **Control** contain raw materials of human origin. The raw materials are from human blood that has been tested and found to be negative for HBVs antigen, HCV antibody, HIV antibody, HIV RNA, HCV RNA, and syphilis. However, handle them with care as they are potentially infectious. No known test method can offer complete assurance that products derived from human sources will not transmit infectious agents.

##### 8.3 Limitations

- Do not use reagent bottles for purposes other than this test.
- Do not separate reagent bottles that connect to **NESCAUTO Cp Auto NS-Prime**. Do not use them in combination with other reagent bottles even if they are the same lot.
- Do not damage or stain the bar codes on bottle labels.
- Do not replenish or mix reagents. Also, do not mix reagents of different bottles even if they have the same lot number.
- Do not use combinations of products with different lot numbers.
- Use of materials other than those listed on "5.0 ADDITIONAL REQUIRED EQUIPMENT" may produce incorrect results.
- Do not recycle bottles. They may be infectious.
- Do not use glassware** for this test. Calprotectin could be adsorbed to the glass, resulting in incorrect results.
- In the case of diarrhea, collection of feces may be inadequate and correct measurement values may not be obtained.
- Use only to assess calprotectin content in human feces. Do not use other specimens.
- Use **NESCAUTO Cp Auto Calibrator** as a calibrator. Other calibrators may produce incorrect results.
- Patients who are taking NSAIDs regularly may have an artificially elevated fecal calprotectin concentration.<sup>4)</sup>

#### 9.0 PRIMARY SAMPLE COLLECTION, HANDLING, AND STORAGE

Use human fecal specimens for measurement.

- To prepare the sample, thoroughly scrape the surface of the fecal specimen with the collector stick from **Specimen Collection Container A**.
- Specimens whose concentrations exceed the upper limit of the calibration curve should be diluted with **FIT NS-Prime Specimen Diluent** and retested.

##### 9.1 Specimen collection

Sample feces onto the notched tip of the collector stick by scraping several parts of the fecal specimen. Place the stick into the collector body once only and fasten the stick tightly in place. See the instruction manual for **Specimen Collection Container A**.

##### 9.2 Sample storage

The collected sample should be refrigerated at a temperature of between 2–8°C until the sample can be tested. Once the sample is received for testing, store between 2–8°C.

##### 9.3 Calprotectin stability after feces sampling

After collection of the sample, it should be measured within 7 days in room temperature conditions (below 25°C). In order to assess human calprotectin stability, four different feces samples collected with the buffer of **Specimen Collection Container A** were examined. After storage for 32 days at –40°C and 8 days at 4, 25, 37 and 50°C, residual ratios of human calprotectin concentration were as follows.

Table 1: Stability of calprotectin in **Specimen Collection Container A**

	–40°C for 32 days	4°C for 8 days	25°C for 8 days	37°C for 8 days	50°C for 8 days
120 $\mu$ g/g	109%	110%	99%	84%	81%
275 $\mu$ g/g	100%	94%	97%	87%	78%
483 $\mu$ g/g	104%	101%	93%	91%	79%
670 $\mu$ g/g	103%	95%	99%	84%	81%

Note: This result is provided for reference only. Calprotectin stability varies depending on the specimen.

#### 10.0 EXAMINATION PROCEDURE

##### 10.1 Preparation of specimens

- Sample feces using **Specimen Collection Container A**. See the instruction manual for the container.
- Shake the container sufficiently to dissolve the feces from the notched tip of the stick.
- Leave the container at room temperature for 1 hour. When measuring on the next day, refrigerate at 2–8°C.

Note: Depending on the specimen, dissolution of calprotectin from feces may be insufficient at 1 hour.

##### 10.2 Assay

Assay procedures for the **NESCAUTO Cp Auto NS-Prime** are established on **Discrete Clinical Chemistry Analyzer NS-Prime**. The analyzer measures calprotectin concentration following the reaction sequence.

Sample	10 $\mu$ L
Reagent 1	140 $\mu$ L
Reagent 2	50 $\mu$ L

Mix the reaction mixture and incubate at 37°C.

Formula to calculate the change of Absorbance: (Abs) is (Am1–As1)–(Am2–As2).

Am1 : Abs of main wavelength at measurement point 1

As1 : Abs of sub wavelength at measurement point 1

Am2 : Abs of main wavelength at measurement point 2

As2 : Abs of sub wavelength at measurement point 2

main wavelength: 540 nm

sub wavelength: 660 nm

measurement point 1: 0.0 minute

measurement point 2: 6.8 minutes

##### 10.3 Calibration curve

**NESCAUTO Cp Auto Calibrator** is used to construct a calibration curve according to the instrument manual. A new calibration must be performed for each new reagent lot. Otherwise, calibration should be performed every 4 weeks.

#### 11.0 CONTROL PROCEDURE

It is recommended to use **NESCAUTO Cp Auto Control** as quality control materials. The values obtained for the quality control materials should not fall repeatedly outside the acceptable ranges. If these control values fall repeatedly outside of the established control ranges, then proper instrument performance should be verified, or recalibration should be performed.

#### 12.0 CALCULATION OF EXAMINATION RESULTS

Specimen values are calculated by **Discrete Clinical Chemistry Analyzer NS-Prime**.

#### 13.0 INTERPRETATION OF RESULTS

Samples with a calprotectin concentration less than 50  $\mu$ g/g are regarded as negative IBD.

For IBD monitoring, samples with a calprotectin concentration more than 118  $\mu$ g/g are regarded as active IBD.<sup>5)</sup>

However, we recommend each laboratory to establish its own reference concentration range. Various cut-off values have been used in the literatures.<sup>5-10)</sup>

## 14.0 PERFORMANCE CHARACTERISTICS

### 14.1 Analytical performance characteristics

#### 14.1.1 Precision

Claimed Precision:  $CV \leq 15\%$

Repeatability and within-laboratory precision were determined according to the CLSI approved guideline EP5-A3. Six samples were tested in duplicate according to assay procedure in a period of over 20 days with two lot reagents. One run per day was performed. The values of repeatability and within-laboratory precision were 1.0–1.9% and 2.1–3.4%, respectively.

#### 14.1.2 Interference

Some studies have been conducted to determine the levels of interference according to the CLSI approved guideline EP7-A2. The levels of interference due to interfering substances were in the range of 90–110% at the following concentrations in feces (The concentration unit definition was substance mass per gram feces):

Erythrocytes (0.188 g/g)  
 Urine (0.5 g/g)  
 White blood cells (10 mg/g)  
 Hemoglobin (0.54 mg/g)  
 Adalimumab (0.32 mg/g)  
 Azathioprine (1 mg/g)  
 Mesalazine (40 mg/g)  
 Magnesium oxide (20 mg/g)  
 Ciprofloxacin hydrochloride hydrate (6 mg/g)  
 Vancomycin hydrochloride (30 mg/g)  
 Aspirin (3 mg/g)  
 Loxoprofen sodium hydrate (1.8 mg/g)  
 Mefenamic acid (15 mg/g)

Erythrocyte concentrations higher than 0.188 g/g may affect the measured value. Avoid blood contamination as much as possible when collecting samples.

### 14.2 Diagnostic performance characteristics

#### 14.2.1 Analytical sensitivity

The difference in the amounts of change in absorbance between calprotectin concentrations of 0 and 100  $\mu\text{g/g}$  is 0.04–0.14.

#### 14.2.2 Analytical accuracy

Analytical accuracy is within  $100 \pm 20\%$  of known concentration in the measurement of a control specimen of known concentration.

### 14.3 Measuring interval

#### 14.3.1 Assay range

Upper limit: 1200  $\mu\text{g/g}$  (CLSI EP6A)

Lower limit:

LOQ  $\leq 40 \mu\text{g/g}$  (CLSI EP17-A2)

Detection limit  $\leq 35 \mu\text{g/g}$  (CLSI EP17-A2)

### 14.4 Correlation

It had been reported that there is a strong correlation between results from NESCAUTO Cp Auto and the reference method.<sup>5)</sup> A linear regression analysis was performed and the following results were obtained:

x : NESCAUTO Cp Auto  
 y : Reference method (ELISA)  
 $r = 0.98$   
 $y = 0.94x + 6.9$   
 $n = 309$

## 15.0 BIOLOGICAL REFERENCE INTERVALS

Table 2 : Statistically result reference interval study

Population	healthy adults
n	119
Range ( $\mu\text{g/g}$ )	4–457
Mean ( $\mu\text{g/g}$ )	69.4
0–95 <sup>th</sup> percentile ( $\mu\text{g/g}$ )	182
Median ( $\mu\text{g/g}$ )	43

This study was conducted with apparently healthy adults (in-house data).

Note: This value is indicative only and may differ from other published values because of differences in methods and in the population being studied.

## 16.0 LIMITATION OF THE EXAMINATION PROCEDURE















- Cuts, abrasions, and other skin lesions should be properly protected with an appropriate waterproof dressing, because blood contamination may influence measurements.
- Appearance changes, such as cloudiness and aggregation, in any of the reagents indicate the possibility of deterioration. Call your local distributor for advice.
- As gold particles may precipitate during storage, mix the Reagent 1 and Reagent 2 combined bottles thoroughly before each use. Mix it by slowly inverting to avoid creating bubbles. If bubbles form, the instrument sensor cannot detect the liquid surface correctly and misunderstands as if the reagent volume increased (reagent from another bottle was added). As the result, the bottle will be deleted from the system. Pipette out any bubbles down to the level of the liquid.
- As with all assays, the results of this test can be influenced by compounds present in some patients' specimens.
- For diagnostic determination, the results obtained from this assay should always be used in combination with a clinical examination, patient medical history, and other findings.

- Procedural directions must be followed exactly, because any modification of the procedure may change the results. Read the instrument instruction manual and use it according to the described usage and operating environment.
- Use of reagents, disposables, or spare parts other than those supplied by the authorized distributor may produce incorrect results.
- Store the reagents according to the storage methods. Do not use them after the expiration date.
- Use fresh feces.
- Take necessary precautions to reduce infection risk when collecting feces.
- This test should not be used to analyze specimens taken from a patient who is menstruating or who has hemorrhoids.

## 17.0 LITERATURE REFERENCES

- Røseth AG. et al. "Assessment of the neutrophil dominating protein calprotectin in feces. A methodologic study". Scand J Gastroenterol. 1992; 27(9): 793–8.
- Lasson A. et al. "The intra-individual variability of faecal calprotectin: a prospective study in patients with active ulcerative colitis". J Crohns Colitis. 2015; 9(1): 26–32.
- Okuyama Y. et al. "A novel sol particle immunoassay for fecal calprotectin in inflammatory bowel disease patients". Clin Chim Acta. 2016; 456: 1–6.
- Tibble JA. et al. "High prevalence of NSAID enteropathy as shown by a simple faecal test". Gut. 1999; 45(3): 362–6.
- Inoue K. et al. "Usefulness of a novel and rapid assay system for fecal calprotectin in pediatric patients with inflammatory bowel diseases". J Gastroenterol Hepatol. 2014; 29(7): 1406–12.
- Lobatón T. et al. "A new rapid test for fecal calprotectin predicts endoscopic remission and postoperative recurrence in Crohn's disease". J Crohns Colitis. 2013; 7(12): e641–51.
- Wright EK. et al. "Measurement of fecal calprotectin improves monitoring and detection of recurrence of Crohn's disease after surgery". Gastroenterology. 2015; 148(5): 938–47.
- Naismith GD. et al. "A prospective evaluation of the predictive value of faecal calprotectin in quiescent Crohn's disease". J Crohns Colitis. 2014; 8(9): 1022–9.
- Lasson A. et al. "Pharmacological intervention based on fecal calprotectin levels in patients with ulcerative colitis at high risk of a relapse: A prospective, randomized, controlled study". United European Gastroenterol J. 2015; 3(1): 72–9.
- Bressler B. et al. "Clinicians' guide to the use of fecal calprotectin to identify and monitor disease activity in inflammatory bowel disease". Can J Gastroenterol Hepatol. 2015; 29(7): 369–72.

## 18.0 SYMBOLS USED IN PRODUCT INSERTS AND ON LABELS

Symbols	Meanings of the symbols
	CE marking
	Use-by date
	Batch code
	Catalogue number
	Manufacturer
	Authorized representative in the European Community (Authorized European representative)
	Number of tests
	<i>In vitro</i> diagnostic medical device ( <i>In vitro</i> diagnostic)
	Temperature limit (for store)
	Consult instructions for use
	Caution: Products containing hazardous substances
	Mixing of substances.
	Reagent 1
	Reagent 2

## 19.0 DATE OF ISSUE OR REVISION

December 1, 2020